


Get better with ActivaNail™

Excellent for small fractures.

BETTER FIXATION

- For fixation of bone fractures, osteotomies, arthrodesis and osteochondral fractures in the presence of appropriate immobilization.
- Improved fixation stability with flake fractures.
- Self-Locking SL™ and grooved surface design enable nails to stay in place with improved rotational stability.

FDA cleared and CE approved



Example of an indication:
OCD (Osteochondritis
Dissecans) with ActivaNail™

EASY TO USE

- High strength properties offer easy insertion and safe medical use.
- Aseptic use.

EXCELLENT PRODUCT PROPERTIES

- Well controlled degradation.
- No stress shielding.
- No need for removal operation.
- Bioabsorption reduces risks of long term complication.
- Implant is supplied Gamma sterilized.
- Manufactured in Finland according to ISO 13485:2003 certified quality management system.

Product codes and instruments

ActivaNail™ Conical

B-ANCO-1510	ActivaNail™ Conical 1.5 mm x 10 mm
B-ANCO-1515	ActivaNail™ Conical 1.5 mm x 15 mm
B-ANCO-1520	ActivaNail™ Conical 1.5 mm x 20 mm
B-ANCO-1525	ActivaNail™ Conical 1.5 mm x 25 mm
B-ANCO-1530	ActivaNail™ Conical 1.5 mm x 30 mm
B-ANCO-2015	ActivaNail™ Conical 2.0 mm x 15 mm
B-ANCO-2020	ActivaNail™ Conical 2.0 mm x 20 mm
B-ANCO-2025	ActivaNail™ Conical 2.0 mm x 25 mm
B-ANCO-2030	ActivaNail™ Conical 2.0 mm x 30 mm



Instruments

B-IP-4000	Applicator handle, arthroscopic, for pins and nails
B-IP-4004	Applicator housing, arthroscopic, for pins and nails
B-IP-3212	Applicator tip, arthroscopic, for 1.5/2.0 mm nails
B-IP-4001	Applicator drill guide, arthroscopic, for 1.5/2.0 mm nails
B-IP-1504	K-Wire Ø1.6 mm, length 260 mm, arthroscopic, for 1.5 mm pins/nails
B-IP-1505	K-Wire Ø1.6 mm, length 260 mm, cone tip, arthroscopic, for 1.5 mm pins/nails
B-IP-3200	Applicator Ø3.2 mm, reusable, for 3.2 mm pins and for 1.5/2.0/2.7/3.2 mm nails
B-IP-2004	K-Wire Ø2.0 mm, length 260 mm, arthroscopic, for 2.0 mm pins/nails
B-IP-1503	Drill bit Ø1.5 mm, for 1.5 mm pins/nails, quick coupling
B-IP-2003	Drill bit Ø2.0 mm, for 2.0 mm pins/nails, quick coupling



Vertrieb:

BlueRock Medical
Possenhofener Str. 5
D-82340 Feldafing

☎ +49 (0) 8157 998 1238
☎ +49 (0) 8157 998 1242
✉ info@bluerockmedical.de

Manufacturer:

bioretect

Bioretect Ltd.

Hermiankatu 22
FI-33720 Tampere, Finland