


Get better with **ActivaScrew™ Cannulated**

Easy insertion with guide wire.




RESULTS THAT MEET REQUIREMENTS IN SURGERIES

- ActivaScrew™ Cannulated implants have been used in more than 15,000 operations all around the world since 2007. The clinical track-record related to implant usage in surgeries is excellent.

RELIABLE AND EASY FIXATION

- No need for removal operation.
- For fixation of bone fractures, osteotomies, arthrodesis, bone grafts and osteochondral fractures in the presence of appropriate immobilization.
- The functionality of the ActivaScrew™ Cannulated begins during surgery, actively maintaining fixation during and after implantation, which is conducive to proper ossification.
- The special characteristic of the Bioretec screws is Auto-Compression™ to reduce the risk of unstable fixation. Auto-Compression™ tightens the screw lengthwise contraction of the screw.



Example of an indication:
Fixation of Tibial Plateau
with ActivaScrew™ Cannulated 4.0 mm.

EXCELLENT PRODUCT PROPERTIES

- High strength and stiffness facilitate stable fixation.
- Compression on the fracture line.
- Optimal strength retention during healing and resorption in tissue within approximately two years.
- No stress shielding.
- Bioabsorption reduces risks of long term complication.
- Packaging, instrument and implant design enables aseptic and easy use in surgical procedure.
- Designed and manufactured in Finland according to .ISO 13485:2003 certified quality management system.

Product codes and instruments


Wide size availability of ActivaScrew™ Cannulated

B-ASC-3520	ActivaScrew™ Cannulated 3.5 x 20 mm	B-ALC-3524	ActivaScrew™ Cannulated 3.5 x 24 mm
B-ASC-3524	ActivaScrew™ Cannulated 3.5 x 24 mm	B-ALC-3526	ActivaScrew™ Cannulated 3.5 x 26 mm
B-ASC-3528	ActivaScrew™ Cannulated 3.5 x 28 mm	B-ALC-3528	ActivaScrew™ Cannulated 3.5 x 28 mm
B-ASC-3532	ActivaScrew™ Cannulated 3.5 x 32 mm	B-ALC-3530	ActivaScrew™ Cannulated 3.5 x 30 mm
B-ASC-3536	ActivaScrew™ Cannulated 3.5 x 36 mm	B-ALC-3535	ActivaScrew™ Cannulated 3.5 x 35 mm
B-ASC-3540	ActivaScrew™ Cannulated 3.5 x 40 mm	B-ALC-3540	ActivaScrew™ Cannulated 3.5 x 40 mm
B-ASC-4040	ActivaScrew™ Cannulated 4.0 x 40 mm	B-ALC-3545	ActivaScrew™ Cannulated 3.5 x 45 mm
B-ASC-4045	ActivaScrew™ Cannulated 4.0 x 45 mm	B-ALC-4035	ActivaScrew™ Cannulated 4.0 x 35 mm
B-ASC-4050	ActivaScrew™ Cannulated 4.0 x 50 mm	B-ALC-4040	ActivaScrew™ Cannulated 4.0 x 40 mm
B-ASC-4055	ActivaScrew™ Cannulated 4.0 x 55 mm	B-ALC-4045	ActivaScrew™ Cannulated 4.0 x 45 mm
B-ASC-4060	ActivaScrew™ Cannulated 4.0 x 60 mm	B-ALC-4050	ActivaScrew™ Cannulated 4.0 x 50 mm
B-ASC-4065	ActivaScrew™ Cannulated 4.0 x 65 mm	B-ALC-4055	ActivaScrew™ Cannulated 4.0 x 55 mm
B-ASC-4070	ActivaScrew™ Cannulated 4.0 x 70 mm	B-ALC-4060	ActivaScrew™ Cannulated 4.0 x 60 mm
B-ASC-4080	ActivaScrew™ Cannulated 4.0 x 80 mm	B-ALC-4065	ActivaScrew™ Cannulated 4.0 x 65 mm
B-ASC-4090	ActivaScrew™ Cannulated 4.0 x 90 mm	B-ALC-4070	ActivaScrew™ Cannulated 4.0 x 70 mm
B-ASC-4540	ActivaScrew™ Cannulated 4.5 x 40 mm	B-ALC-4080	ActivaScrew™ Cannulated 4.0 x 80 mm
B-ASC-4550	ActivaScrew™ Cannulated 4.5 x 50 mm	B-ALC-4090	ActivaScrew™ Cannulated 4.0 x 90 mm
B-ASC-4560	ActivaScrew™ Cannulated 4.5 x 60 mm	B-ALC-4540	ActivaScrew™ Cannulated 4.5 x 40 mm
B-ASC-4570	ActivaScrew™ Cannulated 4.5 x 70 mm	B-ALC-4545	ActivaScrew™ Cannulated 4.5 x 45 mm
B-ASC-4580	ActivaScrew™ Cannulated 4.5 x 80 mm	B-ALC-4550	ActivaScrew™ Cannulated 4.5 x 50 mm
B-ASC-4590	ActivaScrew™ Cannulated 4.5 x 90 mm	B-ALC-4555	ActivaScrew™ Cannulated 4.5 x 55 mm
		B-ALC-4560	ActivaScrew™ Cannulated 4.5 x 60 mm
		B-ALC-4565	ActivaScrew™ Cannulated 4.5 x 65 mm
		B-ALC-4570	ActivaScrew™ Cannulated 4.5 x 70 mm
		B-ALC-4580	ActivaScrew™ Cannulated 4.5 x 80 mm
		B-ALC-4590	ActivaScrew™ Cannulated 4.5 x 90 mm



Instruments

B-ISC-1000 ActivaScrew™ Cannulated Instrument Set
Individual instruments are also available.



B-ISC-1000

ActivaScrew™ Cannulated
Patents: FI 124190. EP, US, FI Patent Pending Version 0.0



Vertrieb:

BlueRock Medical
Possenhofener Str. 5
D-82340 Feldafing

+49 (0) 8157 998 1238
+49 (0) 8157 998 1242
info@bluerockmedical.de

Manufacturer:

bioretec

Bioretec Ltd.

Hermiankatu 22
FI-33720 Tampere, Finland

www.bluerockmedical.de